

PTO/SB/97 (08-00)

Approved for use through 10/31/2002. OMB 0651-0031
U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it contains a valid OMB control number.

TO: E. THERKORN

Tel. 703 308 0362, fax 703 872 9310

Certificate of Transmission under 37 CFR 1.8

I hereby certify that this correspondence is being facsimile transmitted to the
United States Patent and Trademark Officeon 7/18/02
Date

Signature

Ashok K. Shukla

Typed or printed name of person signing Certificate

Note: Each paper must have its own certificate of transmission, or this certificate must identify
each submitted paper.Total no. of Pages 3 incl. ThisTitle: Incision-Based FILTRATION/SEPARATION
PIPETTE TIP

Inventors: Ashok K. SHUKLA et al.

Serial No. 09/591,009

Att. with: 1723

Examiner: E. THERKORN

Filing date: 06/09/00

Official

GROUP 1700

JUN 8 1 2002

FAX RECEIVED

Burden Hour Statement: This form is estimated to take 0.03 hours to complete. Time will vary depending upon the needs of the individual case. Any comments on the amount of time required to complete this form should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, Washington, DC 20231. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Assistant Commissioner for Patents, Washington, DC 20231.

#12/6/02
7/23/02
Official

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Re: Patent Application of: Ashok K. Shukla, et.al.

Serial No. 09/591,009

Art Unit: 1723

Filing Date: 06/09/00

Examiner: THERKORN, E

Title: INCISION-BASED FILTRATION/SEPARATION PIPETTE TIP

Honorable Commissioner of Patents and Trademarks
Patent and Trade Mark Office
Washington DC 20231

Date: July 18, 2002

In response to communication mailed on 04/18/02.

Examiner has rejected our explanation based on comparisons with Valaskovic patent (US 6190,559). But the two devices are entirely different.

- Valaskovic patent has a column where the end of the column from where the slurry is drawn into the column is much bigger (several times the diameter of the particles) than the size of the packing material. Whereas, in our device, one end of the tube has a very fine slit or perforation, much smaller than the size of the packing particle so there is no scope for it to escape the tube even under pressure during sample preparation.
- Valaskovic patent uses capillary action to fill the tube whereas in our device, we load it from the wide open end.
- Valskovic patent uses the method of slow evaporation to pack the column once the slurry has been drawn into the column. In our case, the column is loaded from the wide-open end and the solvent escapes through the perforation or slit much smaller than the size of the particle without any fear of losing the material even under pressure or vacuum application.

FAX RECEIVED
JUN 8 1 2002
GROUP 1723


- We are ready to add in our claim 1 that columns are packed using pressure or vacuum.
- We have talked to several scientists who use our product and they are also of the opinion that this product is very unique and new in the market and no such device is available in the market. If required we can get letters from these scientists.

No new matter has been added and a withdrawal of the rejection is requested.

If you have any further questions, please feel free to call me. (Tel 410 997 0301)

Thanking you in advance for all your cooperation.

Respectfully submitted,


Ashok K. Shukla
10316 Kingsway Court
Ellicott City, MD 21042
Tel: 410-465-2212
Fax: 410-997-0772.

Official

FAX RECEIVED
JUN 8 1 2002
GROUP 1700

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☒ FADED TEXT OR DRAWING
- ☒ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☐ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.